


Natus Newborn Care Supplies

For Caregivers
of Premies and
Full-Term Newborns



Carter ~ born at 28 weeks



Estee & Belen ~ born at 35 Weeks



Olivia ~ born at 24 Weeks





Natus Pediatrics

Bringing You More Products for the NICU and Newborn Nursery

Dear Caregivers,

Thank you for bringing our products into your NICUs and nurseries.

The Natus Newborn Care Family of Product combines with our Hawaii Medical brand of great infant care products to bring you the quality supplies you can count on when caring for your infants. We are proud to offer these innovative products that are clinically superior, better for babies and can help save hospitals money.

Our goal is to ensure that you have the resources you need to improve the health and maximize the comfort of the infants and children in your care. That's why the Natus Pediatrics team focuses exclusively on products for you – caregivers in the NICU, PICU and newborn nursery.

We welcome your feedback, your product ideas and suggestions. Please contact us at www.natus.com. With your support, we can continue to make a difference for babies and caregivers worldwide.

Sincerely,

Natus Pediatrics



NeatNick™

Sweeping Action Heel Lancet

Makes a Neater Nick, Less Painful & Costs Less, Too



NEATNICK 2.0 HAS ALL THE IMPROVEMENTS YOU'VE BEEN WAITING FOR – INCLUDING AN EASY-TO-REMOVE PULL RING AND EASY-TO-SEE ARROW



TWIST-TIP VIALS MEAN YOU CAN ADMINISTER SUCROSE JUST A DROP AT A TIME BEFORE PAINFUL PROCEDURES – LIKE HEEL STICKS

NeatNick is a less painful lancet that offers great blood flow and helps reduce recollection rates.

- Faster blade speed means less pain for babies
- Smooth firing – does not “jump”
- Clean-edged incision
- Controlled depth comparable to other brands
- Full-Term size (Mango Orange) for babies weighing more than 1,800-2,000 grams (4-5lbs)
- Premie size (Lime Green) for babies weighing less than 1,800-2,000 grams (4-5lbs)

ITEM #	DESCRIPTION	QUANTITY
1030083	Premie Size – Lime Green Heel Lancets	Box of 100
1030085	Full-Term Size – Mango Orange Heel Lancets	Box of 100

STERILE BLADE, SINGLE-USE DISPOSABLE SHARP. LATEX-FREE.



The NeatNick – TootSweet combo boxes make it easy and convenient to offer pain management for heel sticks. Available together, one box of TootSweet 24% sucrose solution and one NeatNick heel lancet. For more TootSweet product info, see page 3.

ITEM #	DESCRIPTION	QUANTITY
103008PT1-100	1030083 NeatNick Premie 104002B-100 TootSweet 1ml vial	Box of 100 each Box of 100 each
103008FT2-100	1030085 NeatNick Full-Term 10400214-100 TootSweet 2ml vial	Box of 100 each Box of 100 each

STERILE BLADE, SINGLE-USE DISPOSABLE SHARP. LATEX-FREE.

SugarPlum™

Glucose Lancet

When You Only Need a Drop of Blood



1.0MM YELLOW FOR MICROPREEMIES TO FULL-TERM BABIES. (g26)

1.8MM GRAY FOR FULL-TERM BABIES AND ADULTS. (g23)

SugarPlum is the same lancet style used for adult glucose tests – but sized for babies. Avoid unnecessary incisions in your babies' heels – when just a drop of blood is needed use SugarPlum!

ITEM #	DESCRIPTION	QUANTITY
1040043	For Micropreemies to Full-Term Babies Un-bagged, 1.0mm - 26G depth (yellow)	Box of 100
1040041	For Full-Term Babies and Adults Un-bagged, 1.8mm - 23G depth (gray)	Box of 100

STERILE, SINGLE-USE DISPOSABLE “SHARP” USED TO PUNCTURE AN INFANT’S HEEL FOR COLLECTING A BLOOD SAMPLE. LATEX-FREE.

NeatNick 2.0™
Better Than Ever!

GumDrop™

Pacifier

Finally, a Pacifier That Fits Your Baby's Needs... Without Trimming



- Single-body construction – won't come apart
- Soft medical grade silicone covers entire surface
- Meets Consumer Product Safety Commission standards for pacifiers
- The lowest profile pacifier – accommodates all CPAP and nasal cannulas
- Approved for hospital and home use
- No need to cut-down or modify your pacifiers – cut-down pacifiers violate federal law
- Vanilla scent available – some babies and moms find it soothing

Why GumDrop Stays In Babies' Mouths Longer

- Lighter weight than other silicone pacifiers
- Low profile and curved to fit baby's face
- Slightly textured nipple
- Nipple modeled after the most commonly used baby bottle nipple and pacifier shapes

Pacifier Use Can Help Reduce Incidence Of SIDS

The American Academy of Pediatrics guidelines, published November 2005, recommend offering a pacifier at nap-time and bed-time throughout the first year of life.¹

LIGHTER WEIGHT THAN OTHER SILICONE PACIFIERS – GUMDROP STAYS IN BABIES' MOUTHS LONGER – NURSES AND PARENTS APPRECIATE THAT



PREEMIE SIZE
GREEN



FULL-TERM SIZE
ORANGE



FULL-TERM SIZE
PINK



FULL-TERM SIZE
BLUE

ITEM #	DESCRIPTION	QUANTITY
1040011	Full-Term, Mango Orange	Box of 100
1040012	Full-Term, Mango Orange, Vanilla Scented	Box of 100
1040013	Preemie, Lime Green	Box of 100
1040014	Preemie, Lime Green, Vanilla Scented	Box of 100
10400125	Full-Term, Blue	Box of 100
10400132	Full-Term, Blue, Vanilla Scented	Box of 100
10400124	Full-Term, Pink	Box of 100
10400131	Full-Term, Pink, Vanilla Scented	Box of 100

PRODUCT SPECIFICATIONS:

Polycarbonate plastic inner guard over molded with medical grade silicone. Phthalate-free, latex-free and BPA-free. Each GumDrop™ pacifier is packaged clean and ready for use in its own plastic bag.

It's Easy for Parents To Buy GumDrops

Parents can buy GumDrop pacifiers directly at 1.800.396.1550 or online at: www.gumdroppacifier.com. Or, sell GumDrop at your hospital gift shop and pharmacy.

QUOTES from PARENTS!



"...I tried (another) pacifier at the hospital but it was too heavy... GumDrops are so easy for her to keep in her mouth. Plus the nipple is shaped like a mother's nipple... a big bonus since she's nursing."
~ NICOLE'S MOM, WHITEFIELD, NH

"My daughter got this pacifier... in the NICU... Not only did it help her figure out nursing right away but it is an amazing soothing thing for her. I love it!"
~ AMY, CEDAR RAPIDS, IA

1. THE CHANGING CONCEPT OF SUDDEN INFANT DEATH SYNDROME
PUBLISHED ONLINE, NOVEMBER 1, 2005 PEDIATRICS VOL. 116 No. 5 NOVEMBER 2005, PP.1245-1255

TootSweet™

24 % Sucrose Solution

Economical & Convenient... In Twist-Tip™ Vials or Cups



TWIST-TIP VIALS – CONVENIENT FOR SINGLE DOSES – LESS WASTE – PURPLE VIALS TO INDICATE ORAL USE

Keep Babies Safe and Save Money with Easy-to-Use TootSweet

TootSweet helps calm and soothe babies in distress and during painful procedures, while protecting against bacterial contamination.

- Baby-appropriate preservatives prevent bacterial growth
- 12ml in 1oz cups
- 0.5ml, 1ml and 2ml Twist-Tip vials
- Vials reduce cost, inventory hassles and waste
- Single dose, single patient use
- Economical and easy to use
- No refrigeration required
- New purple vials – purple is becoming a standard to indicate oral use

Does Your Sucrose Pass the Tests for USP Preservative Safety Effectiveness?

TootSweet's formulation of preservatives has been tested in accordance with USP Preservative Effectiveness Tests (27 NF 22 2004). The test results meet current USP criteria for Antimicrobial Preservative Effectiveness.

The FIRST Sucrose Solution That's Guaranteed Bacteria and Mold-Free

TootSweet is the first and only sucrose solution that we know of to meet these standards. That means with TootSweet you can be sure your newborns and fragile preemies are not exposed to potentially harmful sucrose borne bacteria. TootSweet's preservatives are found in Ora-Sweet™, Gentamicin, Caffeine Citrate, and other products commonly administered to preemies and full-term babies.



THE LARGE TOOTSWEET CUP ALLOWS EASY DIPPING OF OUR GUMDROP PACIFIER – ELIMINATES OVERFLOW

ITEM #	DESCRIPTION	QUANTITY
24% Clear Sucrose Solution		
1040021	1oz cups with 12ml	Box of 40
1040022	1oz cups with 12ml (6 boxes of 40)	Case of 240
1040029-100	Twist-Tip vials – 0.5ml each	Box of 100
10400212-100	Twist-Tip Vials – 0.5ml indiv. bagged	Box of 100
1040027-100	Twist-Tip vials – 1ml each	Box of 100
10400213-100	Twist-Tip vials – 1ml indiv. bagged	Box of 100
1040028-100	Twist-Tip vials – 2ml each	Box of 100
10400214-100	Twist-Tip vials – 2ml indiv. bagged	Box of 100

TOOTSWEET IS TESTED IN ACCORDANCE WITH USP PRESERVATIVE EFFECTIVENESS TEST (27 NF 22 2004). HAWAII MEDICAL PRODUCTS ARE LATEX-FREE.

FINALLY! SUCROSE IN CONVENIENT TWIST-TIP™ VIALS



Clinicians told us they wanted a convenient way to administer sucrose for heel sticks and other painful procedures. The Twist-Tip is perfect when you want a small amount of sucrose, when you're pressed for time or when babies are in transport and a cup would be unstable and messy. The smaller dose size means less waste and the Twist-Tip vials cost a few pennies less than the cups.

JUST TWIST OFF THE TIP

- Administer just a drop at a time.
- Each drop is approximately 40ul (microliters)
- Less waste and no more spills



StarBaby™

Beanbag Positioning Aid — Cleanable & Reusable



STAR SHAPE PROVIDES MANY POSITIONING OPTIONS

Meets a wide variety of positioning needs for your babies – use to assist a facilitative tuck, to position head and neck, hips, extremities or as an extra set of hands

- Asymmetrical arms and legs – place and adjust individually
- Approximate weight: 1lb

ITEM #	DESCRIPTION	QUANTITY
1060061-6	StarBaby beanbag positioning aids 12.75" x 12.5"	Box of 6
1060062-24	Disposable covers	Box of 24

HAWAII MEDICAL PRODUCTS ARE LATEX-FREE.

Not for sale in the USA

Shell-O™

Gel Pillow and Positioning Aid — Affordable, Long-Lasting Gel Support



FOLD TAIL UNDER TO THICKEN PILLOW FOR LARGER BABIES

Helps reduce pressure to tender preemie heads – Shell-O can help minimize head molding and skin breakdown in preemies and immobilized patients

- No need for multiple products or sizes – 7 inch Shell-O adapts as your preemies grow
- Rectangle is a perfect mattress for micropreemies; pillow or hip roll for bigger babies
- Lasts significantly longer than other brands – single patient use only
- Makes a great cold or hot pak for nursing mothers

ITEM #	DESCRIPTION	QUANTITY
102002HM	7" Shell Shape	Box of 12
1020021	Disposable Covers (7")	Box of 50
102007HM	6" x 9" Rectangular Cushion	Box of 12
1020071	Disposable Covers (6" x 9")	Box of 50

HAWAII MEDICAL PRODUCTS ARE LATEX-FREE.

OraSwab™

Foam Tip Applicator

Give your babies the comfort of oral care with the gentle foam tip swab applicator



BABIES ENJOY IT

Foam tip applicator is soft, safe and absorbent for your oral care applicator needs

- Available in two convenient sizes: extra-small 6mm wide and small 9 mm wide
- Intended use is general cleaning and removing bacteria with solution
- Cotton under the foam provides extra absorbency

ITEM #	DESCRIPTION	QUANTITY
006481	Extra-Small 6mm wide and small 9mm wide OraSwabs 1 pair per pouch, 200 pouches per box	400 per box

PRODUCT SPECIFICATIONS:

CUSTOM MEDICAL GRADE FOAM OVER COTTON, NON-STERILE. ORASWABS ARE 20MM LONG AND 155MM HANDLE.

LifeGuard™ Electrodes

With Hydrogel

Designed to Last Even in High Humidity Environments



Designed for tender baby skin

- Two-layers of hydrogel – one layer secures electrode to skin; the second ensures electrode bonding
- Patterned mylar backing – easy to see if backing drops on the floor

ITEM #	DESCRIPTION	QUANTITY
1020013	LifeGuard Electrodes Radio translucent	3 Electrodes per Pouch 100 Pouches per Box

HAWAII MEDICAL PRODUCTS ARE LATEX-FREE.

STRONG, RELIABLE ADHESION, EVEN IN HIGH HUMIDITY INCUBATORS

LilyPad™

Pressure Relieving Memory Foam Mattress — Keeps Babies More Comfortable



MEMORY FOAM DISTRIBUTES WEIGHT EVENLY AND REDUCES CAPILLARY BED COMPRESSION. LILYPAD COMPLIES WITH FEDERAL FLAMABILITY STANDARDS

Supple, supportive and soft – LilyPad helps reduce fidgeting and helps minimize sleep disruption

- Reduces pressure-induced soft tissue damage and prevents pressure sores in immobilized patients
- Cleanable cover is treated with Ultrafresh DM50 and meets anti-bacterial and anti-fungal standards
- Economical: Clean and re-use, doesn't stain
- Helps reduce vibration during transport

VARIETY OF DIFFERENT SIZES

- For incubators, radiant warmers and transport mattresses.
- See our product sheet or call for sizes.

HAWAII MEDICAL PRODUCTS ARE LATEX-FREE.

CookieCap™

CPAP Stabilizing Cap — Keep CPAP Tubes Comfortably & Securely in Place



- Gentle, secure attachment at temple reduces pressure on nares
- Chin strap prevents cap from sliding and stretches to allow nipple feeding and pacifiers
- Stretchable, soft, foam-padded material stays in place
- Invented by the father of twin preemies
- Safe – no pins or rubber bands needed
- Open top for scalp vein access and bilirubin therapy

ITEM #	DESCRIPTION	QUANTITY
For Dual Tube CPAP		
1060021-15	Preemie Size (450 to 1,600 grams; 1 to 3.5 pounds)	Box of 15
1060022-15	Infant Size (1,600 to 4,550 grams; 3.5 to 10 pounds)	Box of 15
Replacement Chin Straps		
1060025-50	Chin Strap for Preemie Size	Box of 50
1060026-50	Chin Strap for Infant Size	Box of 50

HAWAII MEDICAL PRODUCTS ARE LATEX-FREE.

AVAILABLE IN PREEMIE AND INFANT SIZES

SunFish™ & Essentials™

Perforated Temp Probe Covers & Disposable Probes

The First 2-in-1 Size Temp Probe Covers to Fit Premies & Full-Term Babies



SUNFISH – UNIQUE PERFORATED DESIGN MEANS ONE SIZE FITS PREMIES AND FULL-TERM BABIES

No More Cutting Temp Probe Covers

- One cover fits full-term newborns and preemies
- Optional locking feature stops temp probes from slipping
- Helps keep temperature readings accurate and uninterrupted
- Lasts for days, even when used to reposition the probe
- Hydrogel firmly but gently adheres to babies' skin



ESSENTIALS – OHMEDA STYLE

ESSENTIALS – AIR SHIELDS STYLE

ITEM #	DESCRIPTION	QUANTITY
SunFish Temperature Probe Covers		
1030051	SunFish perforated temperature probe covers	25 per Pouch 10 Pouches per Box
1030052	SunFish perforated, temperature probe covers with special locking feature	25 per Pouch 10 Pouches per Box
Essentials Disposable Temperature Probes with SunFish Perforated Covers		
1030032	Disposable Temperature Probes with SunFish perforated hydrogel cover. (COMPARES TO OHMEDA PART No. 208-0697-700)	Box of 40
1030042	Disposable Temperature Probes with SunFish perforated hydrogel cover. (COMPARES TO AIR SHIELDS PART No. 68-209-70)	Box of 40

HAWAII MEDICAL PRODUCTS ARE LATEX-FREE.



Exceptional eye protection for neonatal phototherapy

Biliband eye protectors are designed to block light to the baby's eyes during phototherapy. Their unique Y-shaped design allows them to remain securely in place.

ITEM #	DESCRIPTION	QUANTITY
900644	Micro: Head Circumference 20–28 cm	20 units per pack
900643	Premature: Head Circumference 24–33 cm	20 units per pack
900642	Head Circumference 30–38 cm	20 units per pack

MiniMuffs
NEONATAL NOISE ATTENUATORS

Giving newborns the quiet time they need

Designed specifically for premature babies, MiniMuffs Neonatal Noise Attenuators protect their sensitive ears and help to reduce the disruptive noise levels associated with special care procedures & transport.



ITEM #	QUANTITY / DESCRIPTION
040341	36 individually wrapped pairs of MiniMuffs Neonatal Noise Attenuators + Instructions



Save
the
gonads
www.savethegonads.com
U.S. Patent Pending

Superior x-ray protection for neonates

Save the Gonads shields are specifically designed to protect the reproductive organs of babies receiving x-rays. Their unique heart shape allows placement on both boys & girls.

ITEM #	DESCRIPTION	QUANTITY
900835	Nano-Premature	25 units per pack
900772	Micro-Premature	25 units per pack
900773	Premature	25 units per pack
900774	Newborn	10 units per pack
900845	Sample Pack*	10 units per pack

*Sample Pack contains 3 Nano-Premature, 3 Micro-Premature, 3 Premature and 1 Newborn.

neonatal
OXYGEN DELIVERY SYSTEMS

Clearly...a better solution for oxygen therapy

Utilizing a transparent & aerodynamic design, Natus oxygen hoods provide a superior solution for neonatal oxygen therapy.



The Natus Pediatrics Product Family



ALGO® Newborn Hearing Screeners

Bio-logic ABAer® Hearing Screener

Echo-Screen® Hearing Screener

neoBLUE® LED Phototherapy Systems

Olympic Bili-Lite™ & Bili-Bassinet™
Phototherapy Systems

Billiband® Eye Protector

Olympic Cool-Cap® System

Olympic CFM™ 6000

Brainz BRM3

Olympic Brainz Monitor

Olympic Smart Scales®

Olympic Warmettes®

Olympic Papoose Boards™

Olympic Circumstraint™

Olympic Warm-Lamp™

MiniMuffs® Neonatal Noise Attenuators

Save the Gonads™ X-ray Shields

Neonatal Oxygen Hoods

Neometrics Data Management

Hawaii Medical Product Family



“Our daughter, Ellie, was born not breathing for 14 minutes...we were told she'd most likely die within the night. She was transferred (University of Michigan Mott Hospital) to undergo cooling therapy using your product (**Olympic Cool-Cap System**).

Today she is a happy and completely healthy 2.5 year old, with no developmental delays. We cannot express our gratitude enough.”

~ JENNY & CHRIS Y., THE WOODLANDS, TX



“We liked the **Billiband** better than the goggles our first preemie had to use! The goggles didn't stay on nearly as well as the Billiband.”

~ JASON W., GREENWOOD, IN





WHAT PEOPLE ARE SAYING



“Wonderful!
As a neonatal nurse and a mother of a strictly breastfed baby,
GumDrop really is a little miracle paci.”
~ **SUSANA, MIAMI, FL**

“They (**CookieCap**) are the BEST!
I highly recommend these for the NICU ... great with CPAP,
feeding tubes, and nasal canulas.”
~ **ROBI, CHESTERFIELD, SC**



“Superior blood flow, clean incision, reduced “milking” ...
NeatNick (is) reducing the number of PKU re-collections.”
~ **MOUNT CARMEL REGIONAL MEDICAL CENTER LABORATORY DEPT.,
PITTSBURG, KS**

“Fantastic even for breastfed babies!
We got our first **GumDrop** in the hospital when Noah was born – it fits him
better than the bigger pacifiers and he loves it. He has been breastfed since the
beginning and has never had any nipple confusion ... they are great!”
~ **BRIENNE P., SHAWNEE, OK**



PRODUCTS & SPECIFICATIONS INCLUDED IN THIS CATALOG ARE SUBJECT TO CHANGE WITHOUT NOTICE.



natus
pediatrics



Natus Medical Incorporated

1501 Industrial Road • San Carlos, CA 94070 USA | 1.800.255.3901 (US/CA) • Worldwide: 1.650.802.0400
www.natus.com