

natus[®]
pediatrics



Shell-O[™] Gel Pillow and Positioning Aid

Affordable, Long-Lasting Gel Support.

- Unique shell shape offers positioning flexibility
- Use just the corner for micropreemies and for smaller babies
- Fold tail under to thicken pillow for larger babies
- Also makes a great hot or cold pack for nursing moms



Give every baby the benefit of gel pillows.
Shell-O lasts your baby's entire length of stay!

- Helps minimize head molding and skin breakdown
- Provides uniform support for the soft cranium of the pre-term infant
- Lasts significantly longer than other brands
- Unique shell shape offers positioning flexibility: Fold the tail under to adjust the shape and thickness as your babies grow
- Ideal for micropreemies, preemies and full-term babies
- Also available in 6" x 9" rectangular shape
- Disposable, baby-soft covers available for each size
- Single patient use only

Shell-O adapts as your preemies grow!



Hawaii Medical
Shell-O[™]



With Shell-0™, Every Baby Can Benefit from Gel Support That's Affordable and Long-Lasting.

Gel pillows are an effective way to minimize pressure and provide uniform support for the soft cranium of the pre-term infant. The wonderfully soft surface and supple support of Shell-0 Gel Pillows keeps babies more comfortable and minimizes the re-positioning needed to alleviate pressure points. Shell-0 also can aid in the healing process if sores or ulcers are present. The Shell-0 patented film and proprietary gel formula means they last for your baby's entire length of stay... saving you money while giving all babies the benefits of gel support!

Easy to Use, Find and Store.

With its sea blue color, Shell-0 won't get lost in the bed, so you'll save money on replacement pillows. And its patented, liquid-lock, baby soft film minimizes evaporation – so it lasts for months! Because you can adjust the size and thickness of Shell-0, one size can be used for micro-preemies, preemies, and full-term babies... even as they grow. That means inventory requirements are reduced and ordering and stocking are easier.

PRODUCT SPECIFICATIONS:

OUTER MATERIAL: Medical grade polyolefin film


CONTENTS: water, propylene glycol, methyl paraben, neutralized carbopol, food coloring. Product is latex-free.



SHELL-0 IS ALSO AVAILABLE IN A 6" x 9" RECTANGULAR SHAPE – PERFECT MATTRESS FOR MICROPREEMIES OR HIP ROLL FOR BIGGER BABIES.

LASTS SIGNIFICANTLY LONGER THAN OTHER BRANDS! SINGLE PATIENT USE ONLY.

No Need for Multiple Products or Sizes – 7" Shell-0 Adapts as Your Preemie Grows!



FOR SMALLER BABIES, SLIDE THE GEL INTO THE TAIL TO ADJUST THE THICKNESS.



FOLD THE TAIL UNDER FOR BIGGER BABIES TO PROVIDE A THICKER PILLOW.

Ordering Information

Give your babies the benefit of gel pillows with Shell-0™. For questions, to place an order or for the sales representative in your area, please call us.

| ITEM # | DESCRIPTION | QUANTITY |
|----------|-----------------------------|------------|
| 102002HM | 7" Shell Shape | 12 per box |
| 1020021 | Disposable Covers (7") | 50 per box |
| 102007HM | 6" x 9" Cushion | 12 per box |
| 1020071 | Disposable Covers (6" x 9") | 50 per box |



natus
pediatrics



Natus Medical Incorporated

1501 Industrial Road • San Carlos, CA 94070 USA | 1.800.303.0306 (US/CA) • Worldwide: 1.650.802.0400
www.natus.com