

- Lasts for years
- Asymmetrical arms and legs- place and adjust individually
- Economical – can be reused between patients
- Approximate weight: 1 lb.

StarBaby was conceived by a NICU nurse!



STARBABY GIVES YOUR BABY COMFORT AND GIVES YOU A HELPING HAND.

Introducing StarBaby! – a cleanable, reusable bean-bag that meets a wide variety of positioning needs for your babies. Just wipe down between use – no laundering required!

- Cleanable treated urethane cover has anti-bacterial and anti-fungal properties.
- No laundering required – just wipe with hospital cleaning solution.
- More effective and economical than blanket rolls (just ask your hospital about laundering costs)!
- Star shape gives you more options for creativity to meet each baby's need.
- Soft, quiet fabric is developmentally appropriate; breathable yet waterproof.

With StarBaby, you can easily, gently position your baby and provide containment comfort.

StarBaby's asymmetrical arms and legs can be placed and adjusted individually, used as an extra hand during procedures or for facilitative tucking. StarBaby is partially filled with polypropylene beans that can be shifted into virtually endless positions. Soft, quiet fabric is developmentally appropriate, breathable yet waterproof. Urethane fabric and underliner safely contain the beans and provide a waterproof, cleanable surface. Just wipe down with your hospital cleaner between babies.

Ordering Information

Get flexible, easily cleanable positioning aids for your babies with StarBaby. For questions or to order, please call us.

ITEM#	DESCRIPTION	QUANTITY
Cleanable, reusable beanbag 12.75" x 12.5"		
1060061-6	StarBaby beanbag positioning aids.	6
1060062-24	Disposable covers	24
<small>PRODUCT SPECIFICATIONS: StarBaby is MRI compatible. Polyester knit with urethane coating, virgin polypropylene filling. These products are latex-free. Manufactured in the USA.</small>		

Hawaii Medical
StarBaby™

